

COVID-19

Antigen detection (CE/IVD-marked)



Simple
Fast
Cost effective
Accurate
Comfortable
Frequent



COVID-19 antigen detection

- Ensure safety and health
- Keep companies and productions up and running
- Give people peace of mind in difficult times
- Act responsible and be an example to follow
- Comfortable testing. Low risk of damaging glands
- Easy handling. Robust and reliable test

Results within 15 minutes

Our new Covid-19 antigen detection test kit for screening of corona virus gives you a result within 15 minutes. Onsite. The test is inexpensive and easy to use. It's easy to store and requires no specialist equipment or machinery - making it the perfect choice, if you can't wait for days for the results, and if you need to test often to stay ahead of the virus.

The Covid-19 antigen test is CE/IVD-marked.

If safety is the ambition, frequency is the answer

We believe that frequent testing is the key to keeping COVID-19 in check. Our antigen-based COVID-19 tests are nearly as precise as the common known PCR-tests (coincidence rate: 97,3%), but with frequent testing you level out the difference. Instead of relying on day old tests, you get the results today.

Easy to use

The Covid-19 antigen detection test kit is simple and intuitive. To ensure a good sample and a reliable result, the sample must be taken by a health care professional with a swab inserted into both nostrils of the test person. With this robust test you get a reliable result within 15 minutes compared to a PCR test where you risk waiting for 72 hours until getting a result.

From events to facility services


Whether you're interested in testing many people once (at e.g. a concert or soccer match), or you need to test a smaller or larger number of people frequently (at e.g. a production facility), we have a solution for you.

Workflow in brief




STEP 1

Collect nasal sample by inserting swab 2 cm into nostril.



STEP 2



Rotate swab to elute sample and close drip lid.



STEP 3

Add 3 drops onto detection card. Incubate for 15 minutes and interpret test results.

Interpretation of results



Positive: Band is marked on both (C) and (T) line
Negative: Band is marked on the (C) line. No line on (T).
Invalid: No band visible or a band on the (C) line.

For more info visit www.dna-diagnostic.com or contact Tel +45 87 32 30 50 · info@dna-diagnostic.com

**DNA HUMAN
DIAGNOSTIC**

DNA DIAGNOSTIC A/S · VOLDBJERGVEJ 14 · 8240 RISSKOV · DENMARK

About us DNA Human Diagnostic is 100% owned by DNA Diagnostic A/S. DNA Diagnostic A/S is a Danish biotech company founded in 1992, specialized in the development of test kits for Human Health, Food Safety and Animal Health. Our products are sold and used all over the world. At DNA Diagnostic, we believe that the world deserves fast, cost effective and accurate diagnostics.