

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

## 1. Product Identification

**Product names:**

Enterit 4 Swine

**Catalog No.:**

E4S

**Chemical names, common names:**

Not applicable

**Product use:**

The kit is for in vitro screening of swine feces for bacteria

**Synonyms:**

Not applicable

**Manufacturer's name:**

DNA Diagnostic A/S

**Address:**

Voldbjergvej 14  
DK-8240 Risskov  
Denmark

**Business phone:** +45 8732 3050

**Emergency phone:** Use your local number

## 2. Hazards identification

**Classification of the substance or mixture according to 1272/2008 CLP**

One product component have a constituent that requires it is classified; Pre-lysis buffer.

Skin irritation - category 2 – H315

Serious eye damage – category 1 – H318

Chronic aquatic toxicity – category 3 – H412

**Hazard statements**

H315 – causes skin irritation

H318 - causes serious eye damage

H412 - harmful to aquatic life with long-lasting effects

**Precautionary statements**

P264 Wash skin thoroughly after handling

P280 Wear protective gloves and eye protection

P302 + P352 IF ON SKIN: Wash with plenty of water

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

P305 + P351+ P338 +P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON

CENTER

P332 + P313 If skin irritation occurs: Get medical advice

P362 + P364 Take off contaminated clothing and wash it before reuse

P273 Avoid release to the environment

P501 Dispose of contents/container according to national legislation

## Label Elements 1272/2008 CLP



Signal word: Danger

### Other Hazards

The product does not meet the criteria for PBT or vPvB.

## 3. Composition / Information of Ingredients

The Mastit 4 Kit consists of a large box and 1-4 small boxes.

Large box contains:

Four 96-deep well plates with seal and 150 µL Pre-lysis buffer in each well.

Four 96-well plates

Four 96-cap sheets.

Four bottles with 110 mL Wash Buffer

Small box contains:

Four 96-well plates with seal and 15 µL qPCR reaction mix in each tube.

Two bottles with 12 mL Lysis Buffer.

Two vials with 1.3 mL Additive for Lysis Buffer

The following components have a constituent that are either classified, have occupational exposure limits or biological constituents.

Pre-lysis buffer:

CAS/EC-no.	REACH No	Name	Content%	Classification CLP
9036-19-5/-		polyethylene glycol octylphenyl ether	<20%	Acute toxicity – category 4 – oral – H302 Serious eye damage – category 1 – H318 Chronic aquatic toxicity – category 2 – H412

qPCR reaction mix and Additive for Lysis buffer:

CAS/EC-no.	REACH No	Name	Content%	Classification CLP
------------	----------	------	----------	--------------------

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

56-81-5/200-289-5	-	Glycerol	<10%	-
-------------------	---	----------	------	---

## Additional information

See full text of H-phrases in section 16. Occupational limits are listed in section 8, if these are available.

## Special note regarding Biological Constituents

The Biological Constituents of this product are known to be non-infectious and non-hazardous. When this product is used in conjunction with bodily fluids UNIVERSAL PRECAUTIONS SHOULD BE IMPLEMENTED.

## 4. First aid measures

### Emergency overview

Water is the main component of all solutions.

Occupational exposure presents little or no health hazard, however there are exposure limits for glycerol see section 8.

### First aid measures

Basic hygiene should prevent any problems. If contact with this product leads to reddening, inflammation, or other skin or eye irritation, flush the exposures area with running water.

Remove any contaminated clothing, taking care not to contaminate eyes. If chemical is in eyes, open victim's eyes while under gentle running water. Use sufficient force to open eye lids. Have the victim "roll" eyes. As per standard industrial practice, the recommended minimum flushing time is 15 minutes.

The victim should seek immediate medical attention if an adverse reaction has occurred. A copy of label and SDS should be given to the physician or health professional. If chemical is swallowed, and the victim is having a reaction to the product, CALL PHYSICIAN OR POISON CONTROL CENTER FOR MOST CURRENT INFORMATION. If professional advice is not available, induce vomiting. Never induce vomiting or give diluents (milk or water) to someone who is unconscious, having convulsions, or who cannot swallow.

### Eye:

Flush with water in eyewash for at least 15 minutes, holding eyelids open. Remove contact lenses, clean before re-use. Obtain medical attention if symptoms develop (redness, itching, etc.)

### Skin:

Wash with plenty of soap and water for at least 15 minutes. Remove contaminated clothing, which may be washed for re-use. Obtain medical attention if symptoms develop (redness, itching, etc.).

### Inhalation:

Remove person to fresh air. Obtain medical attention unless effects are mild and temporary.

### Ingestion:

Give plenty of water, if conscious. If vomiting occurs naturally, wash out mouth. Be prepared to induce vomiting upon a physician's advice. Obtain medical attention if symptoms develop.

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

---

**Note to Physician:**

Treat symptomatically. No specific antidote.

## 5. Fire Fighting Measures

**Fire extinguishing media**

This Product is not combustible. If this product is in the proximity of a fire, use fire extinguishing materials which are appropriate the materials which are involved in the fire, eg. powder, foam, CO<sub>2</sub> or water fog.

**Unsuitable extinguishing media**

None known.

**Hazardous Combustion products**

Includes carbon dioxide, carbon monoxide, and dense smoke.

**Advice for firefighters**

Standard turnout gear ("bunker gear"). Positive pressure self-contained breathing apparatus. Prevent firefighting water from entering the environment.

## 6. Accidental release measures

**Personal Protective Equipment**

Area ventilation is generally adequate. Chemical-resistant gloves should be worn. Wear standard safety goggles or glasses. All watertight materials and rubber materials offer adequate protection. Avoid contact with eyes and skin.

**Spill cleanup**

Absorb spill. Common absorbent materials should be effective. Deposit in appropriate containers for removal and disposal. Wash area with plenty of water.

**Environmental precautions**

Avoid discharge to lakes, streams, sewers, etc. Accidental releases may be subject to special reporting requirements and other regulatory mandates.

See also sections 8 and 13

## 7. Handling and storage

Storage of some materials is regulated by federal, state, and/or local laws.

**Precautions for safe handling**

See section 8 on exposure controls /personal protection for information on personal protection. Avoid contact with eyes and skin. The product should be used under well-ventilated conditions. Keep closed or covered when not in use.

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

## Storage temperature

Large box: Store at room temperature  
Small box: Frozen, at minus 20 degrees Celsius

## Storage pressure

Ambient.

## Handling procedures

Keep closed or covered when not in use.

## Storage Procedures

Suitable for most general chemical storage areas. Protect from direct sunlight.

## Specific end uses

This product should only be used for applications described in section 1.

## 8. Exposure controls, personal protection

### Control parameters

#### Occupational exposure limits EH40 (Great Britain)

Contained in product constituents qPCR mix and additive for Lysis buffer .

CAS-no:	Name	Limits	Comments
56-81-5	Glycerol, mist	10 mg/m <sup>3</sup> /8 hour TWA	-

DNEL/PNEC: -

### Exposure controls

#### Appropriate technical measures:

Well- functioning area ventilation is generally adequate. No respiratory protection will be needed under normal operating conditions.

#### General information/hygiene measures:

Smoking, consumption of food or liquid, and storage of tobacco, food or liquid, are not allowed in the workroom. Whenever you take a break in using this product and when you have finished using it, all exposed areas of the body must be washed. Always wash hands, forearms and face. Make sure that eyewash and emergency showers are clearly marked. Avoid contact with eyes and skin.

### Personal Protective Equipments

Only CE marked personal protection equipment should be used. Chemical-resistant gloves should be worn (EN 374).

Wear standard safety goggles or glasses (EN 166). No respiratory protection is needed.

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

---

Wear appropriate protective clothing.

## 9. Physical and chemical properties

### For all components of the kit

**Odor**

None

**Boiling point**

Not established

**Melting point or range**

Not established

**pH**

neutral

**Vapor pressure**

Not established but similar to water

**Specific gravity**

Not established

**Solubility in water**

Yes

**Vapor pressure, mm Hg at 20 degrees Celsius**

Not established

**Volatiles**

Not established

**Evaporation rate (water=1)**

Not established

**Viscosity**

Not established

**Appearance and color**

Colorless liquid

## 10. Stability and Reactivity

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

---

## **Reactivity**

Stable under normal and recommended conditions.

## **Chemical Stability**

Stable under recommended conditions for use and storage

## **Conditions to avoid**

Any set of conditions, which should be avoided for water. Direct sunlight. Strong oxidizers and reducers.  
Strong acids and bases. Alkalimetals.

## **Materials with which substance is incompatible**

This product is not compatible with substances which are incompatible with water. Otherwise not established.

## **Hazardous Decomposition Products**

None known

## **Hazardous reactions**

Not expected to occur

## **11. Toxicological information**

Only the pre-lysis buffer component of the kit is classified and the information below is for that constituent.

### **Potential Health Effect**

No long-term health effect from exposure to this product is presently known.

### **Specific target organ toxicity – single exposure**

Not classified

### **Acute toxicity**

Not classified

### **Serious eye damage/irritation**

May cause moderate eye irritation

May cause moderate corneal injury

### **Dermal/skin corrosion/irritation**

Prolong skin contact is unlikely to result in harmful amounts.

Typical for this family of materials.

LD50, Rabbit >2000 mg/kg Estimated

### **Inhalation/respiratory**

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

---

Not classified

## **Oral/ingestion**

Not classified

## **Aspiration hazard**

Not classified

## **Specific target organ toxicity – single exposure**

Not classified

## **Specific target organ toxicity – repeated exposure**

Not classified

## **Carcinogenicity**

Not classified

## **Reproductive toxicity**

Not classified

## **Additional information**

May cause mild temporary irritation in eyes

## **12. Ecological Information**

Only the pre-lysis buffer component of the kit is classified and the information below is for that constituent.

### **Toxicity**

Avoid leakage to sewage system or surface water. Constituent contained in product Pre-lysis buffer:

#### Toxicity to fish

For this family of materials:

Material is toxic to aquatic organisms (LC50 between 1 and 10 mg/l in the most sensitive species)

LC50, Pimephales promelas (fathead minnow), static test, 96 hour, 3.5-4.9 mg/l, OECD test Guideline 203 or equivalent

#### Toxicity to aquatic invertebrates

For this family of materials:

EC50, Daphnia magna (Water flea), 48 Hour, 3.1 mg/l, OECD Test Guideline 202 or equivalent

### **Environmental Fate (Degradation, Transformation, and Persistence)**

Biodegradability: For this family of materials: Material is readily biodegradable. Passes OECD test(s) for ready biodegradability.



# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

---

10-day Window: Not applicable  
Biodegradation: > 60 %  
Exposure time: 28 d  
Method: OECD Test Guideline 301F, B or equivalent

## **Bioaccumulative potential**

No data

## **Mobility in soil**

No data

## **Result of PBT and vPvB assessment**

The product does not meet the criteria for PBT or vPvB

## **Other adverse effects**

None known

## **13. Disposal Considerations**

### **Waste treatment methods**

Avoid discharge to drain and surface water Collect spills and waste in closed leak-proof containers for disposal at the local hazardous waste site. This product, when being disposed of in its unused and uncontaminated state should be treated as a hazardous waste according to EC Directive 2008/98/EC. Any disposal practices must be in compliance with all national and provincial laws and any municipal or local by-laws governing hazardous waste. For used, contaminated and residual materials additional evaluations may be required. Do not dump into any sewers, on the ground, or into any body of water.

### **EWC-code**

15 02 02 absorbents, filter materials (including oil filters not otherwise specified), wiping cloths, protective clothing contaminated by dangerous substances.

## **14. Transport information**

### **This product is not subject to transport regulations**

#### **UN Proper shipping name**

Not determined

#### **UN number**

Not determined

#### **Transport Hazard class**

Not determined

#### **Environmental hazards - MP**

# SAFETY DATA SHEET

Enterit 4 Swine  
DNA Diagnostic A/S

Revised: 20.03.2019  
Replaces: N/A

Not determined

**ID number**

Not determined

**Special precautions for the user:** Not applicable

**Transport in bulk according to Annex II of MARPOL 73/78 and IBC code:** Not applicable

## 15. Regulatory information

**European Union.**

Pressure Equipment (Amendment) Regulations 2011. Chemicals (Hazard Information and Packaging for Supply) Regulations 2009. Control of Substances Hazardous to Health Regulations 2002 (as amended). Merchant Shipping (Dangerous Goods and Marine Pollutants) Regulations 1997. Reporting of Injuries, Diseases and Dangerous Occurrences Regulations 1995 (as amended). Personal Protective Equipment Regulations 2002. Personal Protective Equipment at Work Regulations 1992. Hazardous Waste (England and Wales) Regulations 2005(as amended). EC regulation 19s07/2006 (REACH) Directive 2000/532/EC. Seveso directive: 96/82/EC. EU 453/2010.

**Other countries**

Not determined.

## 16: Other information

**Full test of H-phrases as mentioned in section 3 for constituent of prelysis-buffer**

H302 Harmful if swallowed  
H332 Harmful if inhaled  
H315 Causes skin irritation  
H318 Causes serious eye damage  
H319 Causes serious eye irritation  
H401 Toxic to aquatic life  
H412 Harmful to aquatic life with long-lasting effects

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since DNA Diagnostic A/S cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.