

DNA  
DIAGNOSTIC

# Listeria Velox

## Listeria Velox

Can do it faster...

### Why wait for results?

- **Results in only 20 hours**  
(Incl. enrichment)

- **Multiplex detection**

Choose to detect all Listeria spp. and/or  
L. Monocytogenes

- **Cost efficient**

High throughput and fast product release

- **High sensitivity**

Detects 1 CFU in 25 g sample

Listeria Velox

Certification by

NordVal

ongoing

## Introduction

Listeria is widespread in the environment and it can be found frequently in raw foods. It spreads to processed foods through cross contamination, inadequate processing and poor hygiene. Listeria Velox is a fast and reliable method for detection of *Listeria monocytogenes* and/or *Listeria* species.

Listeria Velox is a rapid *Listeria* qPCR tests for the following categories: Raw and processed meat, poultry and fish/seafood, food swabs and processing environmental samples.

## Advantages

Listeria Velox run in 96 well formats, with a kit size of 384 samples. This makes Listeria Velox products easy to use in high-throughput screenings.

The Velox range operates as an open platform and are validated on machines from leading manufacturers. This gives maximum flexibility for implementation.

Listeria Velox provides a time and cost-saving system for *Listeria* spp. detection and/or *L. monocytogenes*.

## Total Solutions

All DNA Diagnostic kits are supplied with all necessary reagents and consumables, including:

- Enrichment broth
- Plastic consumables incl. Enrichment bags
- DNA extraction kits
- Pre-aliquoted qPCR mastermix

## Other kits in this field

- **Salmonella Velox**  
Kit for 96/960 lysis and qPCR reactions for multiplex detection of *Salmonella* spp. and/or *Salmonella* Enteritidis and Typhimurium.

- **Salmonella Velox 4 feces**  
Kit for 384 lysis and qPCR reactions for multiplex detection of *Salmonella* spp. and/or *Salmonella* Enteritidis and Typhimurium in fecal samples and swabs.

- **Salmonella Velox 4 powder**  
Kit for 384 lysis and qPCR reactions for multiplex detection of *Salmonella* spp. and/or *Salmonella* Enteritidis and Typhimurium in powders, environmental samples and swabs.

## Catalogue numbers

Listeria Velox	LV384D, LV384T
----------------	----------------

For more info visit [www.dna-diagnostic.com](http://www.dna-diagnostic.com) or contact Tel +45 87 32 30 50 · [info@dna-diagnostic.com](mailto:info@dna-diagnostic.com)

**DNA  
DIAGNOSTIC**

DNA DIAGNOSTIC A/S · VOLDBJERGVEJ 14 · 8240 RISSKOV · DENMARK

**About us** DNA Diagnostic A/S is a Danish biotech company established in 1992. DNA Diagnostic develops and manufactures qPCR test kits for rapid identification of pathogenic microorganisms. DNA Diagnostic also makes CE-IVD kits for detecting leukemia related translocations. DNA Diagnostic is ISO 13485 certified.