

DNA  
DIAGNOSTIC

# Salmonella Velox Cocoa

Why wait  
for a **result?**

## Salmonella Velox Cocoa

can do it faster ...

- **Fast Sample to Results**  
Results in 14 hours\*
- **Cost savings**  
In production, time and energy
- **Easy to handle**  
Ease-of-use protocol option
- **High sensitivity**  
Detects 1 CFU in 25 g sample
- **Matrices**  
Chocolate, cocoa powders, cocoa liquor,  
milk-derived products and more

Salmonella  
Velox Cocoa

## Protocol: FAST SAMPLE-TO-RESULTS

- 14 to 20 hour sample to results
- 25g samples to 375g samples
- More manipulation steps
- Release product faster
- Reduce storage (liquid tanks)
- Reduce CO2 footprint
- Reduced media volume (1:5 to 1:9)
- Sensitivity: 1 CFU /sample

## Total Solution

Salmonella Velox Cocoa is supplied with most necessary consumables, including:

- Enrichment broth
- Enrichment bags
- Enrichment additive
- Pre lysis kit
- DNA extraction kits
- qPCR kits

## Protocol: EASE-OF-USE

- 20 hour sample to results
- 25g samples to 375g samples
- Very few manipulation steps
- Equipment needs reduced
- Reduce risk of contamination
- Extended enrichment, hands off
- Reduced media volume (1:5 to 1:9)
- Sensitivity: 1 CFU /sample

## Catalogue numbers

SVC96, SVC960

For more info visit [www.dna-diagnostic.com](http://www.dna-diagnostic.com) or  
contact Tel +45 87 32 30 50 · [info@dna-diagnostic.com](mailto:info@dna-diagnostic.com)

# DNA DIAGNOSTIC

DNA DIAGNOSTIC A/S · VOLDBJERGVEJ 14 · 8240 RISSKOV · DENMARK

**About us** DNA Diagnostic A/S is a Danish biotech company established in 1992. DNA Diagnostic develops and manufactures qPCR test kits for rapid identification of pathogenic microorganisms. DNA Diagnostic also makes CE-IVD kits for detecting leukemia related translocations. DNA Diagnostic is ISO 13485 certified.