

DNA
DIAGNOSTIC

Salmonella Velox

Salmonella Velox

Why wait for results?

- Results in only 5.5 - 9 hours

Foods in 5.5 hours (incl. enrichment) and primary production samples in 9 hours (Incl. enrichment)

- Nordval International certified

Beef, pork, poultry, fish & seafood, primary production samples and environmental testing

- Species differentiation

Choose to detect all *S. spp.* and/or *S. Enteritidis* and *S. Typhimurium*

- Cost efficient

High throughput and fast product release

- High sensitivity

Detects 1 CFU in 25 g sample

Salmonella
Velox
product line

Introduction

The Salmonella Velox product line consists of 3 products; Salmonella Velox, Salmonella Velox S. Enteritidis +S. Typhimurium (SE+ST) and Salmonella Velox Cocoa*.

Salmonella Velox and Salmonella SE+ ST are NordVal certified as rapid Salmonella qPCR tests for testing the following categories; Raw and processed meat, poultry and fish/seafood, food swabs, primary production samples and environmental samples.

*Salmonella Velox Cocoa not Nordval certified.

Advantages

The Salmonella Velox product line has both a shorter enrichment time, and less enrichment media, compared to other conventional Salmonella tests.

This makes Salmonella Velox the fastest certified salmonella test on the market.

Salmonella Velox products run in 96 well formats. All steps are optimized for easy handling. This makes Velox products easy to use in high-throughput screenings.

Salmonella Velox products provides time and cost-saving systems for Salmonella spp. detection and/or S. Enteritidis and S. Typhimurium.


Total Solutions

In general the products in the Salmonella Velox product line are supplied with all necessary consumables, including:

- Enrichment broth
- Enrichment bags
- DNA extraction kits
- qPCR kits

Description

SV960SEST



Catalogue numbers

Salmonella Velox	SV96, SV960, SP384, SF384
Salmonella velox cocoa	SVC96, SVC960
Salmonella Velox SE+ST	SV96SEST, SV960SEST, SP384SEST, SF384SEST

Easy and fast workflow

For more info visit www.dna-diagnostic.com or contact Tel +45 87 32 30 50 · info@dna-diagnostic.com

**DNA
DIAGNOSTIC**

DNA DIAGNOSTIC A/S · VOLDBJERGVEJ 14 · 8240 RISSKOV · DENMARK

About us DNA Diagnostic A/S is a Danish biotech company established in 1992. DNA Diagnostic develops and manufactures qPCR test kits for rapid identification of pathogenic microorganisms. DNA Diagnostic also makes CE-IVD kits for detecting leukemia related translocations. DNA Diagnostic is ISO 13485 certified.