




DNA Diagnostic A/S
Voldbjergvej 14
8240 Risskov
Denmark
Phone. +45 8732 3050
Fax +45 8732 3059
www.dna-diagnostic.com

Certificate of Analysis

HemaVision 922N

Name of Product	HemaVision®-9;22N
Catalogue Number	HV02-922N
Kit lot number	10131.Q3.9.22
Expiration Date	2020-05
Storage conditions	-20°C
Content	<ol style="list-style-type: none">1. M6B PCR-I Primer Mix2. M6B PCR-II Primer Mix3. M8F PCR-I Primer Mix4. M8F PCR-II Primer Mix5. cDNA Primer Mix
Quality Control	The ingredients of the HemaVision kit were tested for functionality and contamination in tests using 1) Leukocyte RNA 2) In Vitro transcribed RNA with specific translocation and 3) H ₂ O as template for cDNA synthesis followed by PCR amplification and visualisation of products by agarose gel electrophoresis.
Quality Control approved by	 Lars Berglund QA DNA Diagnostic A/S

