




DNA Diagnostic A/S
Voldbjergvej 16
8240 Risskov
Denmark
Phone. +45 8732 3050
Fax +45 8732 3059
www.dna-diagnostic.com

Certificate of Analysis

HemaVision 1517N

Name of Product	HemaVision®-15;17N
Catalogue Number	HV02-1517N
Kit lot number	10131.Q3.15.17
Expiration Date	2019-05
Storage conditions	-20°C
Content	<ol style="list-style-type: none">1. M4D PCR-I Primer Mix2. M4D PCR-II Primer Mix3. M8C PCR-I Primer Mix4. M8C PCR-II Primer Mix5. cDNA Primer Mix
Quality Control	The ingredients of the HemaVision kit were tested for functionality and contamination in tests using 1) Leukocyte RNA 2) In Vitro transcribed RNA with specific translocation and 3) H ₂ O as template for cDNA synthesis followed by PCR amplification and visualisation of products by agarose gel electrophoresis.
Quality Control approved by	2015-05  Lars Berglund QA DNA Diagnostic A/S