




DNA Diagnostic A/S
Voldbjergvej 16
8240 Risskov
Denmark
Phone. +45 8732 3050
Fax +45 8732 3059
www.dna-diagnostic.com

Certificate of Analysis

HemaVision 119

Name of Product	HemaVision®-1;19
Catalogue Number	HV03-119
Kit lot number	7131.Q2.1.19
Expiration Date	2018-02
Storage conditions	-20°C
Content	1. PCR Primer Mix 2. cDNA Primer Mix
Quality Control	The ingredients of the HemaVision kit were tested for functionality and contamination in tests using 1) Leukocyte RNA 2) In Vitro transcribed RNA with specific translocation and 3) H ₂ O as template for cDNA synthesis followed by PCR amplification and visualisation of products by agarose gel electrophoresis.
Quality Control approved by	2014-02  Lars Berglund QA DNA Diagnostic A/S